



THE POWER OF LIGHT

Lucibel·le
PARIS
Middle East

OVE MINI

THE OVE MINI FROM LUCIBEL·LE
PARIS PUTS THE BENEFITS OF
LIGHT AT THE HEART OF AN
INNOVATIVE, MOBILE CONCEPT
FOR UNRIVALLED BEAUTY.

OVE MINI's technology is based on skin regeneration using LED photobiomodulation. Red light, a source of beauty and well-being, increases the synthesis of collagen, elastin and hyaluronic acid for a plumping, smoothing and toning effect on the face. To make the most of the OVE MINI's features and avoid any damage caused by incorrect use, we recommend that you read the instructions carefully before use.

**The new era of beauty
through light**

3

P5

SAFETY INSTRUCTIONS

P6

CONTRAINDICATIONS

P7

OPERATION, USE & PROTOCOLS

P8

CARE INSTRUCTIONS &
ENVIRONMENTAL PROTECTION

P9

TECHNICAL CHARACTERISTICS

P10

LEGAL WARRANTY

THE POWER OF LIGHT

THIS BOX CONTAINS:

- 1 OVE MINI
- 1 USB-C cable
- 1 pair of glasses
- 1 storage pocket
- 1 instruction manual



SAFETY INSTRUCTIONS

To avoid any risk of electrocution, fire, electric shock or injury, please follow the following safety measures:

5

WARNING : Do not stare at the lit LEDs. It is essential to wear the glasses supplied because of the intense optical radiation.

Check that the mains voltage is within the range indicated on the 5V / 2A mains adapter.

Only use the accessories supplied by Lucibel-le Paris.

Connect the 5V/ 2A mains adapter to a wall socket that complies with current electrical standards.

Do not use the OVE MINI underwater, in a bath or shower.

Never use the OVE MINI with a damaged cord and/or adapter.

Do not use if the OVE MINI does not switch off at the end of the preprogrammed session or if one side of the OVE MINI is broken or removed.

Do not place the OVE MINI on a damp surface or near a bath, shower or sink.

Use the OVE MINI in a dry place and avoid splashing water. If the OVE MINI falls into water : **do not** touch it and disconnect the mains adapter immediately.

Do not use the OVE MINI under a blanket or similar environment. Excessive heat could be generated, causing a fire, electric shock or personal injury.

Do not drop the OVE MINI.

Do not leave the OVE MINI unattended when it is switched on and plugged in.

Do not open or dismantle the OVE MINI. Any work must be carried out by the manufacturer to avoid any danger. If the OVE MINI is loaned, be sure to enclose the instructions.

WARNING : do not attach the light emitters without wearing the glasses supplied.

Because of the intensity of the light, the glasses supplied must be worn throughout the photobiomodulation sessions.

The OVE MINI **must not** be cleaned or maintained by unsupervised children.

CONTRAINDICATIONS

The **OVE MINI** is suitable for all skin types.

However, it is not recommended for use by :

- Pregnant women
- People with telangiectasia and photo-dermatosis
- People with birthmarks
- People suffering from or at risk of epileptic seizures
- Children

We also **strongly advise** against using the OVE MINI if you are taking photosensitising medication or if you are being treated for an illness that causes photosensitivity.

This device **must not be used** by people with reduced physical, sensory or mental capabilities or who lack the experience and knowledge of how to use the device to understand the dangers involved.

Pre-existing skin changes cannot be attributed to exposure to red light from the OVE MINI.

If any unwanted **side effects** are noted, such as itching within 48 hours of the first use, the user should consult a doctor before using the device again.

The use of red light is not contraindicated with anti-cancer treatments. Although red light is contraindicated in the case of active cancer or cancer in progress, it is, on the contrary, a supportive treatment during anti-cancer treatments.

If in doubt, users should seek their doctor's advice.

OPERATION, USE & PROTOCOLS

GETTING STARTED :

Unwind the USB-C cable and connect it to a 5V / 2A USB mains adapter.

Plug the 5V / 2A USB power adapter into a power outlet.

Wear the protective eyewear supplied.

Press the touch button corresponding to the programme you have chosen.

To interrupt the programme, disconnect the USB-C cable.

7

INSTRUCTIONS FOR USE :

Cleanse the skin before use to optimise light absorption. Before the session, avoid applying any cosmetic products other than those recommended by Lucibel-le Paris, as individual sensitivities to certain cosmetic products can occur.

Sessions should be spaced 48 hours apart (e.g. 1 session on Monday then 1 session on Wednesday). This period corresponds to the time taken for the light energy to be 'absorbed' by the cells.

Optimal effects are achieved after an average of 36 sessions (effects may vary depending on the user).

To assess the effectiveness of OVE MINI and measure the progress of results, it is advisable to take a 1st photo of the treated areas before starting the protocol, a 2nd photo after 18 sessions and a 3rd after 36 sessions.

To benefit from the "delayed effect", take a break of 1 to 2 months before resuming sessions. During this break, the skin will continue to improve.

Plug the OVE MINI into a mains socket (for optimum, safe operation, we advise against plugging into a multiple socket).

Fit protective glasses over your face so that your eyes are fully protected during photobiomodulation sessions.

Apply OVE MINI directly to the skin. Choose the face, body or hair programme. At the end of the session, the device stops **automatically**. Repeat the same process to carry out a new session on the area of your choice.

After 3 consecutive sessions, leave the OVE MINI to rest for 5 minutes before starting a new session. At the end of the sessions, disconnect the OVE Mini from the mains.

OPERATION, USE & PROTOCOLS

FACE :

1 session of 6 minutes per area, every other day for 12 consecutive weeks.

3 areas : left cheek / right cheek / forehead.

The Lucibel.le Paris red LED light has been clinically proven to have anti-ageing effects on the face, helping to :

- reduce the depth of crow's feet wrinkles
- improve skin firmness and elasticity
- even out skin tone
- reduce the diameter of skin pores
- reduce the sagging of the oval of the face
- increase dermal density
- smooth skin texture
- regulate sebum excess

BODY :

1 session of 10 minutes per area every other day for 12 consecutive weeks.

The Lucibel.le Paris red LED light has been shown to have an effect on :

- stretch marks
- orange peel skin

HAIR :

1 session of 7 minutes per area every other day for 8 to 12 consecutive weeks.

The Lucibel.le Paris red LED light has been shown to have an effect on :

- hair quality and density
- hair loss
- stimulating the growth of existing hair or hair with an existing bulb

CARE OF OVE MINI :

After each use, clean the OVE MINI using wipes impregnated with a bactericidal solution.

Do not use any detergents, solvents, chemicals or abrasive sponges.

ENVIRONMENTAL PROTECTION :

OVE MINI cannot be treated as household waste: it must be taken to an appropriate collection point for the recycling of electrical and electronic equipment, which helps to preserve natural resources.

TECHNICAL SPECIFICATIONS

9

en

Operating temperature	+ 10°C to 35°C
Service life	50 000 h
Materials	ABS / PMMA / Verre
Weight	430 g
Dimensions	160 x 90 x 11 mm
Glass size	145 x 85 mm
Number of LEDs	10
Input & Output electrical characteristics	5 Vdc / 2A
Wavelength	635 nm
Power	6.5 W

WARRANTY

The Company guarantees the compliance of Lucibel.le Paris Products with the Company's specifications and all state French, European, and applicable technical standards of each country where its products are marketed.

All Products on sale benefit from the legal guarantee of conformity (articles L.217-4-3 et seq. of the French Consumer Code) and the guarantee against hidden defects (articles 1641 and 1644 of the French Civil Code), allowing the Customer to return defective or non-compliant Products.

The Products sold must conform to the contract and be free from defects in conformity at the time of delivery.

The Vendor shall also be liable for any lack of conformity resulting from the packaging, the assembly instructions or the installation when the latter was its responsibility under the contract or was carried out under its responsibility, or when the incorrect installation, carried out by the consumer as provided for in the contract, is due to shortcomings or errors in the installation instructions supplied by the Vendor.

In accordance with Article L. 217-4 of the French Consumer Code, the goods conform to the contract if they meet the following criteria, where applicable :

- It corresponds to the description, type, quantity and quality, in particular as regards functionality, compatibility, interoperability and any other characteristics provided for in the contract;
- It is fit for any special purpose intended by the consumer, made known to the seller at the latest at the time of conclusion of the contract and accepted by the latter;
- It is supplied with all accessories and installation instructions, which must be provided in accordance with the contract;

It is updated in accordance with the contract. In accordance with article L. 217-5 of the French Consumer Code, in addition to the criteria of conformity with the contract, the product is in conformity if it meets the following criteria :

It is fit for the use normally expected of goods of the same type, taking into account, where applicable, any provisions of European Union law and national law as well as all technical standards or, in the absence of such technical standards, the specific codes of conduct applicable to the sector concerned;

- Where appropriate, it shall possess the qualities which the seller has presented to the consumer in the form of a sample or model, prior to the conclusion of the contract.
- Where applicable, the digital elements it contains are supplied in the most recent version available at the time the contract is concluded, unless the parties agree otherwise;
- Where applicable, it is supplied with all the accessories, including packaging, and installation instructions that the consumer may legitimately expect.

11

However, the Company shall not be bound by any of the foregoing public statements if it demonstrates:

1. That it did not know them and was not legitimately in a position to know them;
2. That, at the time the contract was entered into, the public statements had been corrected under conditions comparable to the initial statements;
3. Or that the public statements could not have influenced the decision to purchase;

Any malfunction resulting from a design, material or manufacturing defect shall entitle the Customer to exchange or repair the Product, after examination of the Product recognised as defective by the Company.

This guarantee applies only to the direct benefit of the end consumer (Customer).

Wheels, power supply cable, and external body damage (including normal wear & tear from use) are not included into the warranty.

The warranty does not cover defects proven to be caused by failure to comply with the instructions for use or maintenance indicated by the Company and included in the user guide.

Similarly, the warranty does not cover defects that are proven to be the result of negligence or a cause unrelated to the Product, such as poor electrical installation, poor maintenance, lack of supervision, fire, water damage, sand or liquid ingress, improper storage, abnormal use, use in a damp or outdoor environment or dropping / falling.

The warranty do not apply in the event of the use of parts or accessories not manufactured by the Company, in the event of alteration or removal of the batch number affixed at the factory, in the event of assembly, disassembly, intervention or modification carried out by a third party not approved by the Company.

To benefit from the warranty, the Customer must refrain from carrying out any repairs themselves and must notify the Company at contact@lucibelleparis.ae of the defects attributed to the Product as soon as these defects become apparent, by providing:

- a detailed description of the alleged defect(s);
- the serial number of the Product and the conditions under which it was used;
- proof of purchase;

The Company will process the warranty claim and notify the Customer.

Once the Company acknowledges receipt of the Warranty Claim and sends a notification authorizing a return, the Customer must send the Product to the Return Address: **STI Innovations LLC, Warehouse 2, Marabea Street, Al Quoz 1, Dubai, UAE** (the complete Product, in its original packaging).

Unauthorised returns may be refused by the Company.

The cost of returning the Product shall be borne by the Customer.

In the event of a defect recognised by the Company under the present warranty, the Company will, at its discretion, replace the Product(s) recognised as defective or exchange the material or parts recognised as defective, with Products of at least equivalent performance. In this case, the cost of reshipping the Product back to the Customer will be borne by the Company.

Any other liability or compensation is excluded. If a Product is replaced under this warranty, the replaced Product is itself warranted, but only for the remainder of the original warranty period.

In no event shall the Company be liable for any direct, indirect, special, incidental or consequential damages (including lost profits) arising from use of the Product not in accordance with this user guide or for purposes other than those for which the Product was designed. In accordance with regulations, this warranty is only applicable to end consumer (the Customer).

SUPPORT

I HAVE A PROBLEM. WHO TO CONTACT?

In case of difficulties using the OVE MINI, first check that the instructions described in this guide have been followed. Ensure that the OVE MINI is properly connected to the power outlet / socket.

If the above has been followed & checked, and the malfunction persists, please contact us by email at the following address:

contact@lucibelleparis.ae

You can also reach us via the following number:

+971 56 115 9343

For safety reasons, do not try to disassemble the OVE MINI yourself. If the device has been taken apart or modified, the warranty will no longer be valid.

REACH OUT TO US FOR PERSONALIZED ADVICE

Purchased the OVE Mask directly on the website?
Not yet in contact with one of our lumino-dermatologists?
Lucibel.le Paris has a dedicated team of beauty professionals to advise you about our technology & products.

To get in touch with us, you are welcome to visit our website **www.lucibelleparis.ae** and fill in the Contact Us form at the bottom of the page.

BEAUTY MASK

OVE

Mini

THE POWER OF LIGHT

Lucibel·le

PARIS

Middle East

OF09-S&C, Arabian Sky Business Center, Citibank Building, Umm Hurair 2, Dubai

www.lucibelleparis.ae



NE-0003

STI INNOVATIONS LIGHTING EQUIPMENT TRADING L.L.C - Main License 1257482 Register No. 3130823 DET DUBAI

